




Product Brief

ConnectedHealth™ Nursing Manager Application

ConnectedHealth's *Nursing Manager Application* ("NMA") automates the collection of medical information for patients in a clinical setting.



Nursing Manager Application Brief

- Automates collection of medical information in a clinical environment, using Continua™ and proprietary medical devices.
- Automatically uploads data to ConnectedHealth Portal.
- Provides a robust and secure execution environment
- Supports two-factor authentication using finger print recognition and NFC
- Supports Bluetooth™ Serial Port Profile, and Health Device Profiles specifications
- Supports Continua™ guidelines.

The NMA automates the capture of data from medical devices such as Glucose Meters, Weighing Scales, Blood Pressure Monitors, and Pulse Oximeters via Bluetooth™ technology. The NMA supports both Continua™ compliant, as well as proprietary devices.

The NMA runs on ConnectedHealth™'s Android Medical Tablet, which provides a robust and secure platform to ensure patient security and privacy.

The system supports advanced biometrics and a secure boot environment, providing a mechanism for clinicians to safely and securely capture and store critical patient information.

Using the Nursing Manager Application, clinicians or nurses can:

- Register patient details and capture patient image and biometric data.
- Validate nurse/clinician access using two-factor authentication (fingerprint and NFC card).
- Utilize fingerprint recognition to validate patient identity prior to capturing data, or utilize integrated camera to correlate patient bar-code information.
- Automatically capture data from medical devices, such as Glucose Meters, Weighing Scales, Blood Pressure Monitors within a clinical setting; and transfer the data to the ConnectedHealth™ eHealth Platform. Clinicians can review historical data utilizing a standalone ConnectedHealth™ Dashboard application, or directly through secure webpages on ConnectedHealth™ eHealth Platform.

SINGAPORE
ConnectedHealth
 Unit 03-20/22
 67 Ayer Rajah Crescent
 Singapore 139950
 Tel: +65 6777 9750
 Fax: +65 6778 9005
 Email: sales@connhealth.com
 Web: www.connhealth.com

USA
Sovran Inc.
 3664 Grand Ave, Suite A
 Oakland, CA 94610
 USA
 Tel: +1 510 645 9377
 Fax: +1 510 645 9379
 Email: sales@sovrans.net

EUROPE
ARS Software GmbH
 Starnberger Str. 22
 D-82131 Gauting
 Munchen, Germany
 Tel: +49 (89) 893 4130
 Fax: +49 (89) 8934 1310
 Email: info@ars2000.com




Product Brief

Supported Wireless Medical Devices

Blood Pressure Monitors	A&D Medical UA-767PBT-C** A&D Medical UA-767 PBT Omron BP792IT (US), 708-BT (EU) FORA D40 2-in-1 BPM and Glucose meter**
Weighing Scales	A&D Medical UC-321PBT-C** A&D Medical UC-321 PBT Omron HBF-206IT (US), BF-206BT (EU)
Pulse Oximeters	Nonin Onyx II 9650**
Glucose Meters	Entra MyGlucoHealth MGH-BT1

Please contact ConnectedHealth™ for the latest list of supported devices ** = Continua™ compliant devices

How it all comes together



Part Numbers

Part Number	Description
CH-411-01	Android Nursing Manager (Proprietary Protocols)
CH-421-01	Android Nursing Manager (Continua™ Compliant)
CH-320-01	Connected Health Android Medical Tablet

All trademarks and trade names are those of their respective owners. ConnectedHealth makes no warranties with respect to this document and disclaims any implied warranties of merchantability and fitness for a particular purpose. Information in this document is subject to change without notice. While every attempt is made to ensure accuracy, ConnectedHealth shall not be liable for any errors or omissions that may appear in this document.